## Checklist upkeep of amalgam separator



Week Number \_\_\_\_\_

EVERY WORKDAY	MON	TUES	WEDS	THURS	FRI		SAT	SUN	DEVIATIONS/OTHER
Disinfectant in the suction system and fountain bowl	S	I G	Ν	Α	Т	ι	JR	Е	
Disinfectant in drains with amalgam separators connected to sinks	S	I G	N	Α	т	L	JR	Е	
Empty coarse strainers/ saliva suction filters/par- ticulate traps into a collec- tion receptacle approved for amalgam waste. Then rinse the filters/traps into in drains with amalgam separators connected to sink.	S	I G	N	Α	т	L	JR	E	
Disruption	S	I G	Ν	Α	Т	L	JR	Ε	
ONCE A WEEK WET SUCTION SYSTEM	MON	TUES	WEDS	THURS	FRI		SAT	SUN	DEVIATIONS/OTHER
Visual inspection of the space in which the amal- gam separator and suction system are installed to check for signs of leakage or noise	S	I G	N	Α	т	L	JR	Е	
Check valve/secretion separator by unfastening the tube ahead of the drain and checking that suction is present at the outlet	S	I G	N	Α	т	L	JR	Е	







Financially supported by EU LIFE Programme

